Zuhause TEST

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EN

INTENDED USE

This product **ZuhauseTEST Syphilis** is a qualitative, manual, *in vitro* diagnostic self-test designed to detect *Treponema pallidum* (TP) antibody in human fingertip blood sample to aid diagnosis of Syphilis infection. This assay provides only a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test.

For users 13-18 years of age, the test should be conducted by a parent or legal guardian.

SUMMARY

Treponema pallidum is the causative agent of human syphilis, mainly sexually transmitted. Clinically, it can be manifested as primary syphilis, secondary syphilis, tertiary syphilis, latent syphilis and congenital syphilis (fetal syphilis). The clinical manifestations are extremely complex. In the early stage, it mainly invades the skin and mucous membranes, and in the late stage, it can invade the blood vessels, central nervous system and various organs of the whole body. Clinically, this product is mainly used for the auxiliary diagnosis of *Treponema pallidum* infection. The traditional laboratory diagnostic techniques are isolation, culture and identification of pathogen. The rapid diagnosis of patients is mainly using chemiluminescence and immunochromatography.

TEST PRINCIPLE

The product **ZuhauseTEST Syphilis** adopts the principle of double antigen sandwich. When the sample contains *Treponema pallidum* antibody, the *Treponema pallidum* antibody in the sample reacts with the colloidal gold-labeled antigen (*Treponema pallidum* recombinant antigen 1) on the conjugate pad to form a labeled antigen-antibody complex, and the complex is chromatographed upward through capillary action, is captured by the detection line (T) antigen (*Treponema pallidum* recombinant antigen 2) coated on nitrocellulose membrane, and a red band appears. The complex continues to be chromatographed upward, and the chicken IgY in the colloidal gold marker is captured by the control line (C) antibody (goat anti-chicken IgY antibody) coated on the nitrocellulose membrane, and a red band appears. When the content of the analyte in the sample is lower than the analytical sensitivity, the detection line (T) does not develop color.

SAMPLE TYPE

The sample type is human fingertip blood.

MAIN COMPONENTS

- 1 Test pad
 - 1 Sample diluent
 - 1 Disposable plastic straw
 - 1 Medical waste bag
 - 1 Blood lancet
 - 1 Alcohol cotton pad
 - 1 Instructions for use
- Note 1:** The components in the kits of different batch numbers are not interchangeable.
- Note 2:** Timer and cotton swab need to be prepared by the users.

STORAGE CONDITIONS AND VALIDITY

Storage conditions: The original packaging should be stored in a dry place at 2-30°C protected from light, and do not freeze. Validity point: 24 months.

The test pad should be used within 1 hour under the condition of 2-30°C and humidity <65% after the aluminum foil bag is unpacked. The production date and expiry date of the product can be viewed on the product labelling.

TEST METHOD

Please read the instructions for use carefully before the test, check whether the device is within the validity period, and

check whether the device is missing or damaged. The device should be returned to room temperature (25°C) before testing, and testing should be performed at room temperature. **Note:** A step-by-step guide can be found on the LAST page and explains how to perform the test.

QUALITY CONTROL

A procedural control is included in the test. The appearance of the control line (C) confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

PERFORMANCE CHARACTERISTICS

- Accuracy, sensitivity and specificity:**
Test results of the candidate kit and the comparative kit are summarized in 2x2 table below.

Zuhause-TEST Syphilis	Comparative device results		
	Positive	Negative	Total
Positive	475	0	475
Negative	2	415	417
Total	477	415	892

Diagnostic sensitivity: 99.58% (95%CI: 98.48% ~ 99.88%)
Diagnostic specificity: 100.00% (95%CI: 99.08%~100.00%)
Overall accuracy: 99.78% (95%CI: 99.19%~99.94%)

ZuhauseTEST Syphilis	Professional test results		
	Positive	Negative	Total
Subject self-testing	429	0	429
Negative	0	370	370
Total	429	370	799

Diagnostic sensitivity: 100.00% (95%CI: 99.11%~100.00%)
Diagnostic specificity: 100.00% (95%CI: 99.97%~100.00%)
Overall accuracy: 100.00% (95%CI: 99.92%~100.00%)

- Interfering substance and cross reactivity:**
Interfering substance: The following potentially interfering substances are added to Syphilis negative and positive samples. None of the substances at the concentration tested interfered in the product.

Rheumatoid factor 5x10² IU/mL, Hemoglobin 5 mg/mL, Bilirubin 20 mg/mL, Procaine penicillin G 0.03 U/mL, Benzathine penicillin 0.03 U/mL, Doxycycline 1 mg/mL, Tetracycline hydrochloride 0.5 mg/mL, Aqueous penicillin G 0.03 U/mL, Procaine 0.5 mg/mL, Ceftriaxone 0.5 mg/mL, Tetracycline 0.5 mg/mL, Erythromycin 1 mg/mL, Sulfamethoxazole 0.5 mg/mL, Isoniazid 1 mg/mL, Rifabutin 1 mg/mL, Pyrazinamide 1 mg/mL, Ganciclovir 1 mg/mL, Acyclovir 1 mg/mL, Valciclovir 1 mg/mL, Sulfadiazine 1 mg/mL, Antibiotic 2 mg/mL, Morphine 100 µg/mL, Trimethoprim 1 mg/mL, Rifampicin 1 mg/mL, Ethambutol 0.5 mg/mL, Clarithromycin 0.35 mg/mL, Sodium fosfarnet 1 mg/mL, Famciclovir 1 mg/mL, Pyrimethamine 2 mg/mL, Glucocorticoid 2 mg/mL, Antifungal drug 2 mg/mL, Antipyretic analgesic 1 mg/mL.

Cross reactivity:
ZuhauseTEST Syphilis does not have the cross reactivity with Measles virus, Rubella virus, Varicella-zoster virus, Epstein-Barr virus, Hepatitis A virus, Hepatitis B virus, Hepatitis C virus, HIV-1 virus, *Helicobacter pylori*, *Mycobacterium tuberculosis*, *Leptospira*, *Mycoplasma pneumoniae*, Influenza A virus, Vaccinia virus, and Human papillomavirus.

- Analytical sensitivity (LoD)**
The analytical sensitivity of **ZuhauseTEST Syphilis** is 3 mIU/mL.

- Hook effect**
ZuhauseTEST Syphilis can test the sample with the concentration of 3072 mIU/mL, and there is no Hook effect.

- Precision**
The precision results of the device within the 20 days all meet the requirements.

WARNINGS AND IMPORTANT INFORMATION

- This product is only used for *in vitro* qualitative diagnosis, please use it within the validity period. There is no devices/equipment needed to be used in combination with the device.
- Please read the instructions for use carefully before use, and operate in strict accordance with the kit instructions for use.
- The product is for single-use. If the aluminum foil bag is found to be damaged, please do not use it. All samples and post-use products should be treated as infectious agents. All samples and post-use products should be properly disposed of in accordance with local relevant regulations.
- If the devices are stored in the refrigerator, it is recommended that they should be taken out of the refrigerator before the test and placed at room temperature before use, otherwise the test results will be affected.
- This device contains animal-derived materials and has potential infectious risks, so try to avoid direct contact with the device.
- Since this product is for visual interpretation, in order to ensure accurate results, please do not interpret in dim light.
- There is a desiccant in the aluminum foil bag, which should not be taken orally; the sample diluent should not be taken orally.
- The sample diluent contains a small amount of preservatives, which may cause irritation to skin and eyes. If this solution comes into contact with skin or eyes, wash/rinse with plenty of water. If skin irritation or rash occurs, seek medical advice/attention.
- The summary of safety and performance about the product can be found on EUDAMED.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- Please guide the result within the specified reading time. Too early or too late may produce invalid results.
- Please read the instructions carefully before operation. Too much or too little sample may produce invalid results.

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LIMITATIONS

- This device is for *in vitro* diagnosis only.
- The test results of this device are for reference only and should not be used as the sole basis for clinical diagnosis and treatment. The clinical management of patients should be comprehensively considered in combination with their symptoms/signs, medical history, other laboratory tests and treatment response.
- Limited by the methodology of the testing device, the testing personnel should pay more attention to the negative results and make a comprehensive judgment based on other testing results. It is suggested that other methods can be used to review the negative results in doubt.
- Weakly positive test results should be confirmed with other methodologies.
- This kit is a reagent for qualitative screening of *Treponema pallidum* antibodies, and the specific concentration of *Treponema pallidum* antibodies in samples cannot be determined.
- Not for use in screening donated blood and tissue for syphilis.
- Persons with a past or current non-venereal *Treponema* infection (Yaws, Pinta) may give positive reactions with this test.

REFERENCES

- Fazal H Z, Gandrakota N, Shackelford J, et al. Mpx infection in an AIDS patient with syphilis manifesting with scrotal and penile cellulitis [J]. Journal of Medical Virology, 2023, 95.
- Shinohara K, Furubayashi K, Kojima Y, et al. Clinical perspectives of *Treponema pallidum* subsp. Endemicum infection in adults, particularly men who have sex with men in the Kansai area, Japan: A case series [J]. Journal of Infection and Chemotherapy, 2022, 28(3):444-450.
- Amaya-Guio J, Grillo-Ardila C F, Edith Angel-Müller, et al. Point of care rapid test for diagnosis of syphilis infection in pregnant women [J]. Cochrane Database of Systematic Reviews (online), 2018 (1).

APPROVAL AND MODIFICATION DATE OF THE SPECIFICATION

Updated version	Updated date	Summary of Changes
1.1	2024.01.05	/
1.2	2024.12.02	Add Hook effect, Analytical sensitivity (LoD) and Precision, etc.
1.3	2025.01.22	Add intended use, some Limitations and Precautions, etc.
1.4	2025.12.19	Translations in German, French, and Spanish have been added, along with information on importers and brands.

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EN Positive: Red bands appear in both the detection line (T) and the control line (C). The results show that the sample contains syphilis antibodies.

Negative: There is no red band in the detection line (T), and a red band appears in the control line (C). The results show

that syphilis antibodies are not detected in the sample.

Invalid: No red band appears in the control line (C), regardless of whether there is a red band in the detection line (T), the test pad is judged to be invalid, and it is recommended to perform a retest.

Note: Do not read the results in dim light. It is not recommended to make any medically related decisions without first consulting a healthcare provider.

If a red background appears in the result display window during the valid reading time, affecting the result reading, it is an invalid result. Please use a new kit to retest or collect samples for retesting. If the problem recurs, please contact the manufacturer through your local distributor at info@nanorepro.com.

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
EN Before the test: Wash your hands before using.

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
EN Open the foil bag, remove the reagent and place it on a horizontal, dry surface.

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
EN Wipe your fingertip clean with the alcohol cotton pad.

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
EN Use the blood lancet to puncture the fingertip. Use the straw to draw the fingertip blood. The first drop of blood may contain too much body fluid, which will affect the test result, so please use a cotton swab to wipe off the first drop of blood.

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
EN Immediately put the straw straightly to add 2 drops of fingertip blood sample, and then add 2 drops of sample diluent into the sample port of the test card.

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
EN Start timing after adding the sample diluent. Read the result in 10–30 minutes and the result is invalid after 30 minutes.

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
EN All used consumables, testing devices and other wastes should be placed in medical waste bags and properly disposed of in accordance with local relevant regulations.

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EN After the test: Wash or disinfect your hands.

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Nanjing Synthgene Medical Technology Co., Ltd.
Building B6-2, No. 9, Weidi Road, Xianlin University Town, Xianlin Subdistrict, Qixia District, Nanjing City, 210000, China
Tel: +86 (025) 83696681
Web: www.syngened.com

CE 2797

REF CR00303001HKS

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NanoRepro AG
MedPath GmbH