de Evolsin® Anti Zahnstein

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en Evolsin® Anti tartar

EVOLSIN® ANTI TARTAR supports the prevention and early removal of tartar in the context of domestic use. It is intended exclusively for use on intact oral mucosa.

1. Indications

- Supports the prevention and early removal of tartar in the context of domestic use.
- Intended exclusively for use on intact oral mucosa and intact gums.
- in people who wear orthodontic appliances and who have soft tissues at risk
- The product is intended for self-application by laypeople in the home environment.

2. Contraindications - The product must not be used:

- in people with inflamed, reddened, swollen or bleeding gums;
- in people with open areas, cracks, ulcerations, or other defects of the oral mucosa;

- in people with blood clotting disorders or who take blood-thinning medication;
- in people who wear orthodontic appliances and who have soft tissues at risk of injury;
- in children under the age of 18 without dental supervision.

3. Warnings and precautions

- Only use on uninjured gums/uninjured oral mucosa
- Only guide the product on the tooth surfaces; do not press or rub it against soft tissue
- Do not use on sensitive teeth
- In the event of pain, bleeding or visible damage to the mucous membrane, discontinue use of the product immediately
- Do not use if the product is damaged
- Do not apply excessive force
- If in doubt about the condition of your gums, consult a dentist beforehand
- Not for daily tooth cleaning. Use only when needed.
- Do not use if you have a known allergic reaction to silicon carbide.

4. Directions

1. Clean your hands thoroughly and check the product
2. Carefully guide the rubbing side along the surfaces of the teeth using very little to no pressure. Only perform light, controlled movements – do not press or scratch.
3. Avoid contact with the gums/mucous membranes
4. Check the pressure again each time you reapply
5. Stop the process as soon as the surface appears smooth
6. Rinse your mouth – do not swallow.
7. Clean and dry the product

5. Possible side effects

- Occasional hypersensitivity of the tooth enamel
 - Temporary discomfort if used improperly
- If you experience the following symptoms, **stop using the product** and seek dental advice:
- Bleeding, persistent pain
 - Injuries to the mucous membrane or gums
 - Visible defects on the tooth

6. Cleaning and storage

- Clean with water after each use
- Store in a dry, clean place
- Do not use the EVOLSIN® ANTI TARTAR device if it is visibly damaged
- Keep out of reach of children

7. Service life

- Replace when the abrasive element is worn

8. Note for patients/users

Any serious incidents related to the product must be reported to the manufacturer and the competent authority of the EU Member State in which you live.

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