Zuhause TEST

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
1. WHO laboratory manual for the examination and processing of human semen - 5th ed. World Health Organization.
2. M.A. Coppola¹, K.L. Klotz², K.-a. Kim³, etc. SpermCheck Fertility, an immunodiagnostic home test that detects normozoospermia and severe oligozoospermia.

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
EN Collect the semen sample into the semen collection cup. If you do lose some of the semen, discard the sample. Rinse the cup **ONLY** with tap water. Wait at least 48 hours, but not more than 7 days and collect a fresh sample.

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
EN Insert the semen transfer syringe into the sample dilution buffer and push the plunger gently to add all of the semen to the solution.

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
EN Gently mix the semen sample in the collection cup and let it stand for **20-60 min**. If the collected semen sample has visibly liquefied, perform the test within 3 hours.

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
EN Screw the enclosed cap back onto the solution bottle and GENTLY mix the contents thoroughly by turning the sample dilution buffer upside down (5-10 times). **DO NOT shake!**

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
EN Use semen transfer syringe to gently stir the semen sample 10 times.

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
EN Insert the semen transfer syringe into the semen sample avoiding any solid or sticky material within the semen. Pull the plunger to draw the sample into the semen transfer syringe until it reaches the **0,1 ml** line on the syringe. Avoid getting air bubbles in the semen transfer syringe. If this happens, push the semen back out completely and then raw semen into the semen transfer syringe again.

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EN Twist off the small cap on the tip of the solution bottle cap and squeeze gently to add **exactly 2 drops** of the semen mixture to the sample well (S) of the test. Please note, that there should be no liquid applied to the result windows marked with the letters (T) and (C).

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Read the result 5-10 minutes after adding the sample. Low values may be missed before the 5 minutes have elapsed, and false positive results may occur after more than 10 minutes.


SCHRITT 5 | STEP 5 | ÉTAPE 5 | PASO 5


SCHRITT 6 | STEP 6 | ÉTAPE 6 | PASO 6


SCHRITT 7 | STEP 7 | ÉTAPE 7 | PASO 7

SCHRITT 8 | STEP 8 | ÉTAPE 8 | PASO 8

SCHRITT 9 | STEP 9 | ÉTAPE 9 | PASO 9

Ergebnis | Result | Résultat | Resultado

Positiv | Positive | Positif | Positivo Negativ | Negative | Négatif | Negativo Ungültig | Invalid | Invalide | No válido